For professional and in vitro diagnostic use only.

The H.pylori Antigen Rapid Test Cassette is a lateral flow immunoassay for the qualitative detection of H.pylori antigen in human feces. It provides an aid in the diagnosis of H.pylori infection.

H. pylori is associated with a variety of gastrointestinal diseases included non-ulcer dyspepsia, duodenal and gastric ulcer and active, chronic gastritis. The prevalence of *H. pylori* infection could exceed 90% in patients with signs and symptoms of gastrointestinal diseases. Recent studies indicate an association of *H. pylori* infection with stomach cancer. *H. pylori* colonizing in the gastrointestinal system elicits specific antibody responses which aids in the diagnosis of *H. pylori* infection and in monitoring the prognosis of the treatment of *H. pylori* related diseases. Antibiotics in combination with bismuth compounds have been shown to be effective in treating active *H. pylori* infection. Successful eradication of *H. pylori* is associated with clinical improvement in patients with gastrointestinal diseases providing a further evidence.

The H.pylori Antigen Rapid Test Cassette is a qualitative membrane strip-based immunoassay for the detection of H.pylori antigen in human feces. In this test procedure, anti-H.pylori antibody is immobilized in the test line region of the cassette. After an adequate volume of test specimen is placed in the specimen well, it reacts with anti-H.pylori antibody coated particles that have been applied to the specimen pad. This mixture migrates chromatographically along the length of the test strip and interacts with the immobilized anti-H.pylori antibody. If the specimen contains H.pylori antigen, a colored line will appear in the test line region indicating a positive result. If the specimen does not contain H.pylori antigen, a colored line will not appear in this region indicating a negative result. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

- For professional in vitro diagnostic use only.
- Do not use after expiration date.
- Do not eat, drink or smoke in the area where the specimens and kits are handled.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Follow standard biosafety guidelines for handling and disposal of potential infective material.
- Humidity and temperature can adversely affect results.

The test contains a membrane strip coated with anti-H.pylori antibody on the test line, goat anti-mouse antibody on the control line, and a dye pad which contains particle coupled with anti-H.pylori antibody.

The quantity of tests was printed on the labeling.

- Test cassette
- Buffer
- Package insert



- Specimen collection container
- Timer

- Store as packaged in the sealed pouch at temperature (4-30°C or 40-86°F). The kit is stable within the expiration date printed on the labeling.
- Once open the pouch, the test should be used within one hour. Prolonged exposure to hot and humid environment will cause product deterioration.
- The LOT and the expiration date were printed on the labeling.

- The H.pylori Antigen Rapid Test Cassette can be performed used on human feces.
- Collect sufficient quantity of feces (1-2 mL or 1-2 g) in a clean, dry specimen collection container to obtain maximum antigens (if present). Best results will be obtained if the assays performed within 6 hours after collection.
- Specimen collected may be stored for 3 days at 2-8°C (36-46°F) if not tested within 6 hours. For long term storage, specimens should be kept below -20°C (-4°F).
- Unscrew the cap of the specimen collection tube, then randomly stab the specimen collection applicator into the fecal specimen in at least 5 different sites to collect approximately 50 mg of feces (equivalent to 1/4 of a pea). Do not scoop the fecal specimen.
- Screw on and tighten the cap onto the specimen collection tube, then shake the specimen collection tube vigorously to mix the specimen and the dilution buffer. Leave the tube alone for 2 minutes.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

Allow the test device and specimens to equilibrate to temperature (15-30°C or 59-86°F) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible.
2. Place the test cassette on a clean and level surface.
3. Holding the specimen collection tube upright, carefully take off the tip of collection tube, then break off the tip of collection tube, transfer 3 drops (approximately 100 μL) to the specimen well (S) of the test cassette, then start the timer. See illustration below.
4. Wait for the colored line(s) to appear. Read results at 15 minutes. Do not interpret the result after 20 minutes.

Positive: Two lines appear. One colored line should be in the control region (C), and another apparent colored adjacent should be in the test region (T). This positive result indicates the presence of *H. pylori* antigens.

Negative: One colored line appears in the control region (C). No line appears in the test region (T). This negative result indicates *H. pylori* antigens are not present or present below the detectable level.

Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test cassette. If the problem persists, discontinue using the lot immediately and contact your local distributor.

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

• The H.pylori Antigen Rapid Test Cassette is for in vitro diagnostic use only. The test should be used for the detection of H.pylori antigen in human feces only. Neither the quantitative value nor the rate of increase in H.pylori antigen can be determined by this qualitative test.

- As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended.

Accuracy

A side-by-side comparison was conducted using the H.pylori Antigen Rapid Test and commercially available H.pylori Ag test. 1068 clinical specimens from three Professional Point of Care sites were evaluated with the H.pylori Antigen Rapid Test and the commercial kit. The following results are tabulated from the clinical study:

		Commercial H.pylori Ag Test		Total
		Positive	Negative	
	Positive	467	1	468
	Negative	2	598	600
Total		469	599	1068

The agreement between these two test kits is 99.57% for positive

specimens, and 99.83% for negative specimens. This study demonstrated that the H.pylori Antigen Rapid Test is substantially equivalent to the commercial device.

Cross-Reactivity and Interference

1. Cross-reactivity with following organisms has been studied at 1.0x10⁹ organisms/mL. The following organisms were found negative when tested with the H.pylori Antigen Rapid Test.

<i>Candida albicans</i>	<i>Campylobacter jejuni</i>
<i>Clostridium difficile</i>	<i>Escherichia coli</i>
<i>Escherichia coli</i>	<i>Enterococcus faecalis</i>
<i>Enterobacter aerogenes</i>	<i>Klebsiella pnemoniae</i>
<i>Proteus vulgaris</i>	<i>Proteus mirabilis</i>
<i>Pseudomonas aeruginosa</i>	<i>Staphylococcus aureus</i>
<i>Salmonella choleraesuis</i>	<i>Salmonella typhi</i>
<i>Salmonella typhimurium</i>	<i>Shigella dysenteriae</i>
<i>Shigella flexneri</i>	<i>Shigella boydii</i>

2. Potentially interference endogenous substances including common components, such as lipids, hemoglobin and bilirubin, were spiked at high concentrations into the H.pylori antigen positive and negative specimens and tested, separately. No interference was observed to the test kit.

Analytes	Conc.	Specimens	
		Positive	Negative
Albumin	20 mg/mL	+	-
Bilirubin	10 µg/mL	+	-
Hemoglobin	15 mg/mL	+	-
Glucose	20 mg/mL	+	-
Uric Acid	200 µg/mL	+	-
Lipids	20 mg/mL	+	-

3. Some other common biological analytes were spiked into the H.pylori antigen positive and negative specimens and tested, separately. No significant interference was observed at the levels listed in the table below.

Analytes	Conc.	Specimens	
		Positive	Negative
Acetaminophen	200 µg/mL	+	-
Acetoacetic Acid	200 µg/mL	+	-
Acetylsalicylic Acid	200 µg/mL	+	-
Benzoylecgonine	100 µg/mL	+	-
Caffeine	200 µg/mL	+	-
EDTA	800 µg/mL	+	-
Ethanol	1.0%	+	-
Gentisic Acid	200 µg/mL	+	-
β-Hydroxybutyrate	20,000 µg/mL	+	-
Methanol	10.0%	+	-
Phenothiazine	200 µg/mL	+	-
Phenylpropanolamine	200 µg/mL	+	-
Salicylic Acid	200 µg/mL	+	-

Reproducibility

Reproducibility studies were performed for H.pylori Antigen Rapid Test at three physician office laboratories (POL). Sixty (60) clinical specimens, 20 negative, 20 borderline positive and 20 positive, were used in this study. Each specimen was run in triplicate for three days at each POL. The intra-assay agreements were > 99%. The inter-site agreement was > 99%.

[BIBLIOGRAPHY]

1. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.

2. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.

3. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996, 100: 35S-41S.

4. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996, 91: 1112-1115.

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Soll, AH. Pathogenesis of peptic ulcer and implications for ther apy. New

England J. Med. (1990), 322: 909-16.

Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Dete ction of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastr oenterology. (1987), 82(4): 292-96.

Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996, 100: 35S-41S.

Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of pep tic ulcer in normal individual with Helicobacter pylori infection. Am J Gastr oenterol. 1996, 91: 1112-1115.